BROEN-LAB QUICK CONNECT WITH BUILD-IN FLOWREGULATION FOR 5.0 GASES

**BROEN-LAB QUICK CONNECT** 



# **BROEN-LAB**

### **BROEN-LAB QUICK CONNECT** Table of contents

Table of Contents	2	
Advantages	2	
Product Range	<b>() () 3</b>	
Outlets	🔖 🤹 🏹 4	
	1.4	
Technical Specifications	5	
Accessories	🝓 🖏 🖏 🛴 5	
Mounting Information	6-7	

# BROEN-LAB QUICK CONNECT Advantages

### Advantages for end users:

- Flexible and compact design increases size of "working area".
- · Reduced downtime.
- Innovative and compact design.
- Combined flow valve and Quick Connect.
- Key system protects against undesired mix of media, 11 media keys are available.
- · Clear media indication.
- Possibility of attaching additional equipment: - flow meters, pressure regulators, etc.
- Non-return valve in outlets protects against contamination of equipment.
- Valve is "locked" with no plug inserted.
- Integrated parking position.
- O-ring which seals outlet and Quick Connect is replaceable without disassembling the Quick Connect.

### Advantages for architects & furniture manufacturers:

- · Increased flexibility in laboratory design.
- · Last minute determination of media possible.
- Number of fittings can be reduced due to the flexibility of the Quick Connect.

### Advantages for plumbers:

- Simple and rational installation.
- 6 different methods of installation.
- Last minute change of media possible (Prior to start of use).
- Possibility to install additional Quick Connect in the laboratory, using a model for exposed piping.

# BROEN-LAB QUICK CONNECT Product range



	Description	Compressed Air	Low Vacuum	Oxygen	Nitrogen	Carbon Dioxide	Hydrogen	Argon	Helium
	Label: According to standard EN13792:2002			02				Ar	He
Model	Label: American recommendation	Air	Vac	02	N <sub>2</sub>	CO <sub>2</sub>	H <sub>2</sub>	Ar	He
1	Standard EN13792:2002	070 <u>020</u> 1	070 <u>030</u> 1	070 <u>040</u> 1	070 <u>050</u> 1	070 <u>060</u> 1	070 <u>070</u> 1	070 <u>080</u> 1	070 <u>090</u> 1
1	American recommendation	070 <u>025</u> 1	070 <u>035</u> 1	070 <u>045</u> 1	070 <u>055</u> 1	070 <u>065</u> 1	070 <u>075</u> 1	070 <u>085</u> 1	070 <u>095</u> 1
2	Standard EN13792:2002	070 <u>020</u> 2	070 <u>030</u> 2	070 <u>040</u> 2	070 <u>050</u> 2	070 <u>060</u> 2	070 <u>070</u> 2	070 <u>080</u> 2	070 <u>090</u> 2
2	American recommendation	070 <u>025</u> 2	070 <u>035</u> 2	070 <u>045</u> 2	070 <u>055</u> 2	070 <u>065</u> 2	070 <u>075</u> 2	070 <u>085</u> 2	070 <u>095</u> 2
2	Standard EN13792:2002	070 <u>020</u> 3	070 <u>030</u> 3	070 <u>040</u> 3	070 <u>050</u> 3	070 <u>060</u> 3	070 <u>070</u> 3	070 <u>080</u> 3	070 <u>090</u> 3
3	American recommendation	070 <u>025</u> 3	070 <u>035</u> 3	070 <u>045</u> 3	070 <u>055</u> 3	070 <u>065</u> 3	070 <u>075</u> 3	070 <u>085</u> 3	070 <u>095</u> 3
Δ	Standard EN13792:2002	070 <u>020</u> 4	070 <u>030</u> 4	070 <u>040</u> 4	070 <u>050</u> 4	070 <u>060</u> 4	070 <u>070</u> 4	070 <u>080</u> 4	070 <u>090</u> 4
4	American recommendation	070 <u>025</u> 4	070 <u>035</u> 4	070 <u>045</u> 4	070 <u>055</u> 4	070 <u>065</u> 4	070 <u>075</u> 4	070 <u>085</u> 4	070 <u>095</u> 4

3 additional media options available on request.

# **BROEN-LAB QUICK CONNECT** Outlets

nci. non- o. 07 <u>XXX</u> 0	with serrated hose nozzle ereturn valve) 0 BROEN-LAB CA		-return va	YNPT female lve)	e	(incl.	ght outlet non-return '0 <u>XX</u> 00	a valve)	
	See dimensions on page 7	1/4" NPT fe	male	See dimensio	ons on page 7	1/4" NP	PT male	See dime	nsions on page
		No. <b>07<u>XXX</u></b>				range: Y flow rang	Pressu	2 nin 0-15 l/min re regulator ber designate	
	See dimensions on page 7	V.		See dimensio	ons on page 9	range: z press. ran Pressure r	<b>3</b> ge 0-3 b	<b>4</b> ar 0-10 bar ressure range of 0-	
	See dimensions on page 7	Compressed Air	Low Vacuum	See dimension		range: z press. ran Pressure r 0-100kPa	<b>3</b> ge 0-3 b egulator with a p	<b>4</b> ar 0-10 bar ressure range of 0 ilable on request	
						range: z press. ran Pressure r 0-100kPa	<b>3</b> ge 0-3 b egulator with a p (0-14.5psi) is ava	<b>4</b> ar 0-10 bar ressure range of 0 ilable on request	-1 bar/
Model	Description Label: According to		Vacuum			range: Z press.ran Pressure r 0-100kPa Carbon Dioxide	<b>3</b> ge 0-3 b egulator with a p (0-14.5psi) is ava	4 ar 0-10 bar ressure range of 0 ilable on request Argon	-1 bar/ Helium
Model 🏹	Description Label: According to standard EN13792:2002 Label: American	Air	Vacuum	Oxygen	Nitrogen	range: Z press.ran Pressure r 0-100kPa Carbon Dioxide	ge 0-3 b egulator with a p (0-14.5psi) is ava Hydrogen	4 ar 0-10 bar ressure range of 0- ilable on request Argon	-1 bar/ Helium
Model	Description Label: According to standard EN13792:2002 Label: American recommendation	Air CA Air	Vacuum Vacuum	Oxygen	Nitrogen	range: z press.ran Pressure r 0-100kPa Carbon Dioxide	ge 0-3 b egulator with a p (0-14.5psi) is ava Hydrogen	4 ar 0-10 bar ressure range of 0 ilable on request Argon Argon	-1 bar/ Helium
Model	Description Label: According to standard EN13792:2002 Label: American recommendation Standard EN13792:2002	Air (CA) Air 07 <u>051</u> 00	Vacuum Vac 07 <u>052</u> 00	Oxygen (0) (0) (0) (0) (0) (0) (0) (0)	Nitrogen	range: z press.ran Pressure r 0-100kPa Carbon Dioxide Coo	ge 0-3 b egulator with a p (0-14.5psi) is ava Hydrogen (H) (H2) (J2) (0705600	4 ar 0-10 bar ressure range of 0 ilable on request Argon Argon Ar O7 <u>057</u> 00	-1 bar/ Helium (He) 07 <u>058</u> 0 07 <u>558</u> 0
Model	Description Label: According to standard EN13792:2002 Label: American recommendation Standard EN13792:2002 American recommendation	Air CA Air 0705100 0755100	Vacuum Vac Vac 07 <u>052</u> 00 07 <u>552</u> 00	Oxygen (0) (0) (0) (0) (0) (0) (0) (0)	Nitrogen (N) (N) (N) (N) (N) (N) (N) (N)	range: z press.ran Pressure r 0-100kPa Carbon Dioxide Coo 0705500 0755500	$\frac{3}{ge} \qquad 0-3 b$ equilator with a p (0-14.5psi) is avain Hydrogen (1) (1) (1) (1) (2) (2) (2) (2) (2) (2) (2) (2	4 ar 0-10 bar ressure range of 0 ilable on request Argon (Ar) (Ar) 0705700 0755700	-1 bar/ Helium Heium 070580 075580 070581
Model	Description         Label: According to standard EN13792:2002         Label: American recommendation         Standard EN13792:2002         American recommendation         Standard EN13792:2002	Air (CA) Air 0705100 0755100 0705110	Vacuum Vac Vac 07 <u>052</u> 00 07 <u>552</u> 00	Oxygen (0) (0) (0) (0) (0) (0) (0) (0)	Nitrogen (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (Ne) (N	range: <b>z</b> press.ran <i>Pressure r</i> <i>0-100kPa</i> <b>Carbon</b> Dioxide <b>Coo</b> <b>0</b> <b>0</b> <b>0</b> <b>0</b> <b>0</b> <b>0</b> <b>0</b> <b>0</b>	$ \begin{array}{c} 3 \\ ge & 0-3 \\ egulator with a p \\ (0-14.5psi) is ava \end{array} $ Hydrogen $ \begin{array}{c} Hydrogen \\ \hline H_2 \\ 0705600 \\ 0755600 \\ 0705610 \\ \end{array} $	4           ar         0-10 bar           ressure range of 0- ilable on request         Argon           Argon         Argon           O705700         0755700           0705710         0705710	-1 bar/ Helium Heium 070580 075580 070581 075581
Model	Description         Label: According to standard EN13792:2002         Label: American recommendation         Standard EN13792:2002         American recommendation         Standard EN13792:2002	Air (CA) Air 0705100 0755100 0705110 0755110	Vacuum Vac Vac 0705200 0755200 07055210	Oxygen (0) (0) (0) (0) (0) (0) (0) (0)	Nitrogen Nitrogen No No No No No No No No No No	range: z press.ran Pressure r 0-100kPa Carbon Dioxide CO2 0705500 0755500 0705510 0755510	$ \begin{array}{c} 3\\ ge & 0-3 b\\ egulator with a p\\ (o-14.5psi) is ava \end{array} $ Hydrogen $ \begin{array}{c} Hydrogen\\ \hline H_2\\ 0705600\\ 0755600\\ 0755610\\ 0755610\\ 0755610\\ \end{array} $	4           ar         0-10 bar           ressure range of 0- ilable on request           Argon           (Ar)           (Ar)           (Ar)           (O705700           0705710           0755710	-1 bar/ Helium Heium O7 <u>0580</u> 07 <u>0580</u> 07 <u>0581</u> 07 <u>0581</u>
Model	Description Label: According to standard EN13792:2002 Label: American recommendation Standard EN13792:2002 American recommendation Standard EN13792:2002 American recommendation	Air (CA) (CA) (Air) 0705100 0755100 0755110 0755110 0705110 070513Y	Vacuum Vac Vac 0705200 0755200 07055210	Oxygen () () () () () () () () () ()	Nitrogen (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N.) (N	range: <b>z</b> press.ran <i>Pressure r</i> <i>0-100kPa</i> <b>Carbon</b> Dioxide <b>CO</b> 2 0705500 0755500 0705510 0705510 0707500 0707500	3         ge       0-3 b         egulator with a p         (0-14.5psi) is avail         Hydrogen         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)         (1)	4           ar         0-10 bar           ressure range of 0-           ilable on request           Argon           (Ar)           (Ar)<	-1 bar/ Helium Hei 070580 075580 075581 075581 070780 070583
Model	Description         Label: According to standard EN13792:2002         Label: American recommendation         Standard EN13792:2002         American recommendation         Standard EN13792:2002         American recommendation         Standard EN13792:2002         Standard EN13792:2002         Standard EN13792:2002         Standard EN13792:2002	Air CA Air 0705100 0755100 0705110 0755110 0707100	Vacuum Vac Vac 0705200 0755200 07055210	Oxygen ()) ()) ()) ()) ()) ()) ()) ()	Nitrogen (N) (N) (N) (N) (N) (N) (N) (N)	range: Z press. ran Pressure r 0-100kPa Carbon Dioxide Co 0 0 0 0 0 0 0 0 0 0 0 0 0	$ \begin{array}{c} 3\\ ge & 0-3 b\\ egulator with a p\\ (o-14.5psi) is ava \\ Hydrogen \\ \hline Hydrogen \\ \hline H_2 \\ 0705600 \\ 0755600 \\ 0755600 \\ 0755610 \\ 0707600 \\ \end{array} $	4           ar         0-10 bar           ressure range of 0- ilable on request         Argon           Argon         Argon           O705700         0755700           0705710         0755710           07070700         0755710	-1 bar/ Helium Helium O70580 070580 070581 075581 070780 070780 070583 075583

# **BROEN-LAB QUICK CONNECT** Technical specifications

### **Quick Connect**

Operating pressure:	-1 bar to 10 bar -14.5 psi to 145 psi
Purity:	5.0 (Only obtainable if system is purged when reconnected)
Test pressure:	1.5 x operating pressure
Operating	0°C to 60°C.
temperature:	23°F to 140°F
Standard media	: Inactive gases
	(Air, Vac, $O_2$ , $N_2$ , $CO_2$ , $H_2$ , Ar & He)
	Non-standard media must be confirmed by BROEN-LAB.
Inlet:	1/4" NPT female
Materials in	
media contact:	Stainless Steel, brass & PVDF
Gaskets:	FKM
External parts:	Stainless Steel, chemical nickel plated brass & EPDM rubber

Outlets

Materials: Stainless Steel and chemical nickel plated brass Gaskets: FKM


### **Pressure regulator**

Recommended P1=2xP2+1 (lower pressure [P1] is possible but will result in minimum inlet lower achievable flowrate) pressure [P1]:

# **BROEN-LAB QUICK CONNECT** Flow charts


### QC + 90° outlet - No. 07<u>XXX</u>00

Example flow chart for Quick Connect and 90° outlet (CA)




### Pressure regulator





0-10 bar regulator



# **BROEN-LAB QUICK CONNECT**

# **Accessories**

Description	Characte	Product No.	
BROEN-LAB UniFlex™ elbow fitting	<b>Material:</b> Brass	<b>Connections:</b> Male 1/4″NPT - Female BROEN-LAB UniFlex™	02101270
BROEN-LAB UniFlex™T-fitting	<b>Material:</b> Brass	<b>Connections:</b> Male 1/4″NPT - 2 x Female BROEN-LAB UniFlex™	02101220
BROEN-LAB UniFlex™ straight fitting	Material: Brass	<b>Connections:</b> Male 1/4″NPT - Female BROEN-LAB UniFlex™	02550220
BROEN-LAB UniFlex™ hose	<b>Material:</b> Braiding: Stainless Steel Plastic: Teflon (PTFE)**** Rubber: FKM**** Metal: Brass	Connections: OD10mm-UniFlex <sup>™</sup> OD12mm-UniFlex <sup>™</sup> OD15mm-UniFlex <sup>™</sup> Union nut G1/2-UniFlex <sup>™</sup> UniFlex <sup>™</sup> -UniFlex <sup>™</sup>	93R1020xxxx*** 93R1220xxxx 93R1520xxxx 93R4020xxxx 93R2020xxxx
BROEN-LAB Flexible Stainless Steel pipe	Material: Braiding: Stainless Steel Metal: Stainless Steel	<b>Connections:</b> Male 1/4"NPT-OD10mm OD10mm-OD10mm	9357080xxxx*** 9358080xxxx

For more information about hoses and their technical specifications please require document 99G0005: BROEN-LAB Uniflex - Flexible Pipes & Hoses

Description	Characteristics		Product No. (all media)
Swagelok** straight tube fitting	Material: Stainless Steel Connections: Male 1/4"NPT	Sube diameter**:           3 mm:           4 mm:           6 mm:           8mm:           1/4":           1/8":	96101010 96101020 96101110 96101030
Swagelok** elbow tube fitting	Material: Stainless Steel Connections: Male 1/4"NPT	Tube diameter: 3 mm: 4 mm: 6 mm: 8 mm: 1/4": 1/8":	96101070 96101080 96100430 96101090
Swagelok** hose connector	Material: Stainless Steel Connections: Male 1/4"NPT	Hose ID: 1/4":	96100450

\* Swagelok – TM Swagelok Company.

\*\* Other dimensions available on request.

\*\*\* Standard lengths of hoses (xxxx):

0500 =500mm, 0700 =700mm, 1000 =1000mm, 1200 =1200mm, 1500 =1500mm, 2000 = 2000mm, 2500 = 2500mm, 3000 = 3000mm


\*\*\*\* PTFE and rubber are penetrable materials meaning that gases, vapors, and liquids may diffuse through them. Always check compatibility with gases in application!


# **BROEN QUICK CONNECT** Mounting information

#### Built-in: No. 070XXX1 Bench/suspended mounted: No. 070<u>XXX</u>4 Ø75 [2 15/16 in] 30 [1 3/16 in] 49 [1 15/16 in] 87 [3 13/32 in] Ø72 [2 27/32 in] Mounting hole Ø52 [2 1/16 in] **®** Ð Ø57 [2<sup>'</sup>1/4 in] ۲ A Ø22 [7/8 in] ANSI B.1.20.1 1/4" NPT 52 [2 1/32 in] Ø51 [2 in] 120 [4 23/32 in] 5 [3/16 in] Panel Max 30 [1 3/16 in] ANSI B1.20.1 1/4" NPT Ø25 [1 in] Max. 32 [1 1/4 in] Table/ceiling Wall/panel mounted, concealed piping 070XXX2. MOUNTED WITHOUT A WALL BRACKET Mounting holes (frontside view) 33.5 [1 5/16 in] M3x0.5 034[1:34] Ø52 2:05 Holes should be bored with a help 12 [0.47] 4.5 [0.18] 0.97 51 [ 2.02] of a wall bracket Mounting hole Contour Column MOUNTED WITH A WALL BRACKET Mounting holes (frontside view) 4.5[0.18] 12 [0.47] 28 [1 3/32 in] Ø52 [2.05] 18 [23/32 in] Hole should be 13 [1/2 in] bored with a help of a wall bracket 023 [0.91] 57 [2.26]


### Wall/panel mounted, exposed piping $\textit{070}\underline{XXX}$ 3. Mounted from the frontside

#### MOUNTED WITHOUT A WALL BRACKET







# **BROEN QUICK CONNECT** Mounting information




### Pressure regulator:

No. 07<u>XXX</u>2<u>ZSP</u>







<b>BROEN</b>	QUICK	CONNE	СТ
Notes			


BROEN	QUICK	CONNEC	Γ
Notes			


<b>BROEN</b>	QUICK	CONNE	СТ
Notes			


BROEN-LAB develops, manufactures and sells Laboratory Fittings, Emergency Shower Systems and Eye Wash Systems into a broad spectrum of laboratories and industrial locations; sectors include Pharmaceutical, Food&Beverage and Academia. Our expertise and product quality ensure optimal solutions compliant to all relevant international norms setting new standards in risk mitigation in modern research and development facilities.

BROEN-LAB offer solutions that ensure the functionality, compliance, hygiene, durability and safety are of the highest priority. For over 50 years our products have been integrated into a wide variety of workplaces, laboratories, hospitals and industrial locations, with features that are assessed and designed to be flexible, durable and compliant offering a broad range of options tailored to each location. This inherent design integrity provides confidence in the solution adopted wherever in the world these are applied.

BROEN-LAB is a collaboration partner in the design and layout of your next laboratory and it's integrated safety solutions.

#### **HEAD OFFICE**

#### **BROEN-LAB A/S**

Drejervaenget 2 • DK-5610 Assens Denmark Tel.: +45 6376 6376 • Fax +45 6471 2476 E-mail: lab@broen.com



#### **BROEN-LAB A/S ISO 9001 certification**

In September 1991 BROEN-LAB A/S was certified according to ISO 9001 as one of the first Danish companies. The certification was carried out by Bureau Veritas, Denmark, for our Danish site in Assens. The quality management system of BROEN-LAB A/S now complies with detailed specifications laid down by the internationally acknowledged EN ISO 9001:2015. This certification will further contribute to reputation for quality and reliability of BROEN-LAB A/S.

#### **INTERNATIONAL OFFICES**

BROEN-LAB Ltd Tel. +44 121 522 4515 E-mail: lab@broen-lab.com

#### BROEN-LAB Singapore Pte. Ltd.

25 International Business Park #04-60A German Centre Singapore 609916 Tel. +65 6298 0662 • Fax +65 6298 0468 E-mail: lab@broen-lab.com

#### **BROEN-LAB GmbH**

Tel.: +45 6376 6376 • Fax +45 6471 2476 E-mail: lab@broen-lab.com

#### BROEN-LAB Sverige AB

Tel.: +45 6376 6376 • Fax +45 6471 2476 E-mail: lab@broen-lab.com

#### **BROEN-LAB Inc.**

15 Constitution Drive Suite 122 Bedford • NH 03110 • USA Tel.: +1 603 310 5089 E-mail: lab@broen-lab.com 046 - BROEN-LAB OUICK CONNECT WITH BUILD-IN FLOWREGULATION- 1

www.BROEN-LAB.com

