


BROEN-LAB

Drejervaenget 2, 5610 Assens, Denmark

Phone: +45 6376 6376

E-mail: lab@broen-lab.com

Web : broen-lab.com



FITTINGS FOR FRONT CONTROL VALVES

PREMIUM range

99B0040 - 30.06.2026


TABLE OF CONTENTS

BROEN-LAB UniFlex™				Page 3
Configuration overview				Page 4
Water		Valve Page 5		Valve Page 6
				Outlet Page 7
Special water		Valve Page 10		
				Outlet Page 11
Coolant water		Valve Page 13		
				Outlet Page 14
Burning gas		Valve Page 15		
				Outlet Page 16
Technical 2.0 gases / vacuum / high flow		Valve Page 17		Valve Page 18
				Outlet Page 19
Technical 4.0 gases / Oxygen		Valve Page 20		
				Outlet Page 21
Pure 5.0 gases		Valve Page 22		
				Outlet Page 23
Pure 6.0 gases		Valve Page 24		
				Outlet Page 25
Corrosive gases		Valve Page 26		
				Outlet Page 27
Accessories				Page 28
Mounting instructions - Front Control Valves				Page 30
Mounting instructions - Outlets				Page 31
Mounting instructions - Mounting of outlets into fume hoods with a double wall construction				Page 32
Pipe installation - Water, Special water (SS), Burning gases, Technical gases 2.0, 4.0 and 5.0				Page 33
Pipe installation - Special water (PP)				Page 33

BROEN-LAB UNIFLEX™ PREMIUM

THIS IS BROEN-LAB UNIFLEX™ PREMIUM

BROEN-LAB UniFlex™ line offers a very high quality and durability combined with the complete program, which ads up to low LCC (Life Cycle Costs). Our line takes the standards in the market to a higher level.




THIS IS BROEN-LAB UNIFLEX™

BROEN-LAB UniFlex™ is a unique and patented connection, which allows customers to choose from various connection methods, such as hoses, CU-, PEX or stainless steel pipes, without changing the valve itself. It provides unique flexibility, not only for the installer, but also for the end-user. Due to the flexibility of the system, the quick installation and savings on fittings, this system offers lower TIC (Total Installed Costs).


PRODUCT OVERVIEW

OUTLETS



FRONT CONTROL VALVES

BROEN-LAB UniFlex™ connection fittings and adaptors
for CU, SS, PEX and PA pipes



For more information about BROEN-LAB UniFlex™ flexible pipes and hoses, connection fittings and adaptors, please contact your BROEN-LAB representative or refer to www.BROEN-LAB.com

FRONT CONTROL VALVE

WATER


Specifications
Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792.
Headwork: Ceramic
Opening/closing function: 180°
Temperature range: 0-90°C
Max. working pressure: 10 bar / 147 psi
Test pressure: 1 x working pressure
Inlet/outlet: BROEN-LAB UniFlex™
Weight: 0.7 kg

Materials with media contact
Metals: Brass
Rubber: EPDM
Plastics: POM
Lubricant: Silicone based
Ceramics: Al-oxyd 96%

Ordering information
For standard configuration please refer to the item numbers below.

- Other options:
- a valve for different medium
 - a valve for front side mounting
 - a valve without mounting kit
 - national approvals

For further information, please contact your BROEN-LAB representative.



Available with or without service valve according to product numbers listed below

Configuration options



Service valve

	without service valve	with service valve
	2511 507 10XX	2511 508 10XX

Replace **XX** with one of these media types:



01 • Potable water, cold



02 • Potable water, hot



07 • Non-potable water, cold



08 • Non-potable water, hot

Delivered with mounting accessories

FRONT CONTROL VALVE

WATER

Specifications

Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792.
Cartridge: Ceramic (with flow and temperature regulation or two step opening)
Temperature range: 0-90°C
Max.working pressure: 10 bar / 147 psi
Test pressure: 1 x working pressure
Inlet/outlet: BROEN-LAB UniFlex™
Weight: 0.6 kg

Materials with media contact

Metals: Brass
Rubber: EPDM
Plastics: POM
Lubricant: Silicone based
Ceramics: Al-oxyd 96%

Ordering information

For standard configuration please refer to the item numbers below.

Other options:
- a valve without mounting kit

For further information, please contact your BROEN-LAB representative.



Configuration options



	with flow and temperature regulation cartridge	with two step opening
	2511 529 10XX	2511 531 10XX

Replace **XX** with one of these media types:

 **35** • Potable water, tempered  **37** • Non-potable water, tempered

Delivered with mounting accessories

OUTLET
WATER

Specifications

BROEN-LAB UniFlex™ outlets have an integrated media indication according to EN 13792.
Inlet/outlet: BROEN-LAB UniFlex™
Max. wall thickness: 20mm
Weight: 0.5 kg

Materials with media contact

Metals: Brass
Rubber: EPDM
Lubricant: Silicone based


Ordering information

For standard configuration please refer to the item numbers below.

Other options:



- other media
- flange
- non-return valve
- fixed nozzle
- other standout
- straight rear wall connection
- wall thickness 20-40mm
- double wall mounting
- table mounting
- additional support in case of a thin panel
- national approvals

For further information, please contact your BROEN-LAB representative.





Additional mounting information see page 30-33

With removable nozzle



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
	2511 101 10XX-40	2511 111 10XX-40	2511 121 10XX-40	2511 131 10XX-40

With quick release coupling



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
	2511 102 10XX-40	2511 112 10XX-40	2511 122 10XX-40	2511 132 10XX-40

Replace **XX** with one of these media types:



01 • Potable water, cold



07 • Non-potable water, cold



35 • Potable water, tempered



02 • Potable water, hot



08 • Non-potable water, hot




37 • Non-potable water, tempered

Delivered with 90° Rear Wall Connection (max wall thickness: 20mm)


OUTLET SPOUT
WATER

With swivel spout




	U = 150 mm	U = 200 mm	S = 150 mm	S = 200 mm
	2511 212 10XX-40	2511 213 10XX-40	2511 214 10XX-40	2511 215 10XX-40

Bench mounted outlet with swivel spout



	U = 150 mm	U = 200 mm
	2511 300 10XX-40	2511 301 10XX-40

45° Outlet with fixed or swivel spout



	Fixed spout	Swivel spout
	2511 307 10XX-40	2511 308 10XX-40

Replace **XX** with one of these media types:

-  **01** • Potable water, cold
-  **07** • Non-potable water, cold
-  **35** • Potable water, tempered
-  **02** • Potable water, hot
-  **08** • Non-potable water, hot
-  **37** • Non-potable water, tempered

Delivered with 90° Rear Wall Connection (max wall thickness: 20mm)

OUTLET SPOUT

WATER

Specifications

Inlet: BROEN-LAB UniFlex™
Hose nozzle: Metal, Removable
Temperature range: 0-65°C
Weight: 0.7 kg

Materials with media contact

Metals: Brass
Rubber: EPDM
Lubricant: Silicone based


Ordering information

For standard configuration please refer to the item numbers below.

Other options:

- non-return valve
- removable nozzle
- other inlet connections (length, type)

For further information, please contact your BROEN-LAB representative.



Ordering information

	2539 039 10XX-52

Replace **XX** with one of these media types:



01 • Potable water, cold



07 • Non-potable water, cold



35 • Potable water, tempered



02 • Potable water, hot



08 • Non-potable water, hot



37 • Non-potable water, tempered

Delivered with mounting accessories

FRONT CONTROL VALVE

SPECIAL WATER

Specifications


Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792.
Headwork: Diaphragm type
Opening/closing function: 1 1/2 x 360°
Max. working pressure: 10 bar / 147 psi
Test pressure: 1 x working pressure
Weight: 0.7 kg

PP valve

Temperature range: 0-65°C
Inlet/outlet: Push couplings for OD10mm pipes

SS Valve

Temperature range: 0-90°C
Inlet/outlet: BROEN-LAB UniFlex™



Available in PP or Stainless Steel according to product numbers listed below

Materials with media contact

PP valve

Diaphragm: PTFE
Plastics: PP

SS Valve

Diaphragm: PTFE
Metals: SS (AISI 316)

Ordering information

For standard configuration please refer to the item numbers below.

Other options:

- a valve without mounting kit
- a valve for front side mounting

For further information, please contact your BROEN-LAB representative.

Configuration options



	PP valve	SS valve
	2511 533 10XX	2511 535 10XX

Replace **XX** with one of these media types:



03 • Distilled water



36 • Deionized water, cold

Delivered with mounting accessories

OUTLET
SPECIAL WATER

Specifications

BROEN-LAB UniFlex™ outlets have an integrated media indication according to EN 13792.
Inlet: Push couplings for OD10mm pipes
Max. wall thickness: 20mm
Weight: 0.5 kg

Materials with media contact

Rubber: EPDM
Plastics: PP


Ordering information

For standard configuration please refer to the item numbers below.


Other options:

- flange
- wall thickness 20-40mm
- swivel spout (U- or S-shape)
- additional support in case of a thin panel

For further information, please contact your BROEN-LAB representative.



Configuration options



	S = 53 mm	S = 75 mm	S = 90 mm
	2511 116 10XX-40	2511 126 10XX-40	2511 136 10XX-40

Replace **XX** with one of these media types:



03 • Distilled water



36 • Deionized water, cold

Delivered with 90° Rear Wall Connection (max wall thickness: 20mm)

FRONT CONTROL VALVE SPECIAL WATER

Specifications

BROEN-LAB UniFlex™ outlets have an integrated media indication according to EN 13792.
Inlet: BROEN-LAB UniFlex™
Max. wall thickness: 20mm
Weight: 0.5 kg

Materials with media contact

Rubber: EPDM
Metals: Stainless Steel (AISI 316)

Ordering information


For standard configuration please refer to the item numbers below

- Other options:
- flange
 - straight rear wall connection
 - fixed nozzle with external G3/4" thread
 - wall thickness 20-40mm
 - additional support in case of a thin panel

For further information, please contact your BROEN-LAB representative.



Configuration options



	S = 53 mm	S = 90 mm
	2511 117 10XX-40	2511 137 10XX-40

Replace **XX** with one of these media types:

- 
03 • Distilled water
- 
36 • Deionized water, cold


Delivered with 90° Rear Wall Connection (max wall thickness: 20mm)

FRONT CONTROL VALVE

COOLANT WATER

Specifications

Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792.
Headwork: Rubber
Opening/closing function: 2 x 360°
Temperature range: 0-90°C
Max. working pressure: 10 bar / 147 psi
Test pressure: 1 x working pressure
Inlet/outlet: BROEN-LAB UniFlex™
Weight: 0.7 kg



Materials with media contact

Metals: Brass
Rubber: EPDM
Lubricant: Silicone based

Available with or without service valve according to product numbers listed below

Ordering information

For standard configuration please refer to the item numbers below.

- Other options:
- a valve for front side mounting
 - a valve without mounting kit
 - a valve with simple hand operation of flow and return of chilled water

For further information, please contact your BROEN-LAB representative.

Configuration options



	without service valve	with service valve
	2511 503 10XX	2511 504 10XX

Replace **XX** with one of these media types:

 **05** • Coolant water, feed  **06** • Coolant water, return

Delivered with mounting accessories

OUTLET
COOLANT WATER

Specifications

BROEN-LAB UniFlex™ outlets have an integrated media indication according to EN 13792.
Inlet/outlet: BROEN-LAB UniFlex™
Max. wall thickness: 20mm
Weight: 0.5 kg
Quick release coupling: double shut-off

Materials with media contact

Metals: Brass
Rubber: EPDM
Lubricant: Silicone based


Ordering information

For standard configuration please refer to the item numbers below.



Other options:

- flange
- other standout
- straight rear wall connection
- wall thickness 20-40mm
- double wall mounting
- additional support in case of a thin panel

For further information, please contact your BROEN-LAB representative.





With removable nozzle
and integrated
non-return valve



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
	2511 104 10XX-40	2511 114 10XX-40	2511 124 10XX-40	2511 134 10XX-40

With quick release
coupling and non-return
valve (return only)



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
	2511 105 10XX-40	2511 115 10XX-40	2511 125 10XX-40	2511 135 10XX-40

Replace **XX** with one of these media types:



05 • Coolant water, feed



06 • Coolant water, return

Delivered with 90° Rear Wall Connection (max wall thickness: 20mm)

FRONT CONTROL VALVE

BURNING GAS

Specifications

Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792.
Headwork: ceramic, press/turn
Opening/closing function: 90°
Max.working pressure: 7 bar / 101,5 psi
Test pressure: 1.5 x working pressure
Inlet/outlet: BROEN-LAB UniFlex™
Weight: 0.7 kg

Materials with media contact

Metals: Brass
Rubber: NBR
Lubricant: Mineral oil based



Available with or without service valve according to product numbers listed below

Ordering information

For standard configuration please refer to the item numbers below.

- Other options:
- a valve for different medium
 - a valve for front side mounting
 - a valve without mounting kit

For further information, please contact your BROEN-LAB representative.

Configuration options



Service valve

	without service valve	with service valve
	2511 519 10XX	2511 520 10XX

Replace **XX** with one of these media types:



10 • Natural Gas



12 • Liquefied Gases



14 • Butane



16 • Propane



18 • Acetylene



20 • Hydrogen



46 • Methane

Delivered with mounting accessories

OUTLET BURNING GAS

Specifications

BROEN-LAB UniFlex™ outlets have an integrated media indication according to EN 13792.

Inlet/outlet: BROEN-LAB UniFlex™

Max. wall thickness: 20mm

Weight: 0.5 kg

Materials with media contact

Metals: Brass

Rubber: NBR

Lubricant: Mineral oil based


Ordering information

For standard configuration please refer to the item numbers below.



Other options:

- other media
- flange
- other standout
- straight rear wall connection
- wall thickness 20-40mm
- double wall mounting
- table mounting
- additional support in case of a thin panel

For further information, please contact your BROEN-LAB representative.





With fixed nozzle



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
	2511 100 10XX-40	2511 110 10XX-40	2511 120 10XX-40	2511 130 10XX-40

With quick release coupling



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
	2511 102 10XX-40	2511 112 10XX-40	2511 122 10XX-40	2511 132 10XX-40

Replace **XX** with one of these media types:



10 • Natural Gas



12 • Liquefied Gases



14 • Butane



16 • Propane



18 • Acetylene



20 • Hydrogen



46 • Methane

Delivered with 90° Rear Wall Connection (max wall thickness: 20mm)

FRONT CONTROL VALVE

TECHNICAL GASES 2.0 & VACUUM/HIGH FLOW

Specifications

Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792.
Inlet/outlet: BROEN-LAB UniFlex™
Weight: 0.7 kg
Test pressure: 1.5 x working pressure

Technical 2.0 gases

Headwork: Needle
Opening/closing function: 3 x 360°
Max. working pressure: 16 bar / 232 psi

Headwork: Micro flow
Opening/closing function: 7.5 x 360°
Max. working pressure: 16 bar / 232 psi

Vacuum / High Flow*

Headwork: High flow
Opening/closing function: 1.5 x 360°
Max. working pressure: -1 bar to 16 bar / -14.5 psi to 232 psi

Materials with media contact


Metals: Brass
Rubber: FKM
Plastics: PVDF
Lubricant: Perfluoropolyether based

Ordering information

For standard configuration please refer to the item numbers below.

- Other options:
- a valve for front side mounting
 - a valve without mounting kit

For further information, please contact your BROEN-LAB representative.



Available with or without service valve according to product numbers listed below

*Please note that valves for vacuum as a standard delivered with a high flow headwork. High flow headwork can also be used for other media where there is a requirement for higher flow.

Configuration options



Service valve

	Without service valve Needle headwork	Without service valve Micro flow headwork	With service valve Needle headwork	With service valve Micro flow headwork
	2511 511 10XX	2511 546 10XX	2511 512 10XX	2511 547 10XX

Replace **XX** with one of these media types:

- 21** • Compressed air
- 23** • Nitrogen
- 24** • Carbon dioxide
- 25** • Argon
- 26** • Helium

	Without service valve High flow headwork	With service valve High flow headwork
	2511 515 10 28	2511 516 10 28

- 28** • Low vacuum (1000 to 1 mbar)

Delivered with mounting accessories

FRONT CONTROL VALVE PRESSURE REGULATOR

Specifications

BROEN-LAB UniFlex™ pressure regulator for panel mounting with a built-in filter (mesh size 200 µm) on the inlet.
Type: Non-relieving
Proof pressure: 15 bar
Operating pressure: Max. 10 bar
Regulating pressure range: 0,5-8.5 bar
Inlet/outlet: BROEN-LAB UniFlex™
Max. panel thickness: 3,5mm
Weight: 0.5 kg

Materials with media contact


Metals: Brass, zinc die-cast
Rubber: NBR, FKM
Plastics: POM
Lubricant: Perfluoropolyether based

Ordering information

For standard configuration please refer to the item numbers below.



Flow rate

Input pressure p1 = 7 bar



Pressure characteristics

Input pressure p1 = 7 bar



Ordering information

	2599 366 1099

21 • Compressed air 23 • Nitrogen 24 • Carbon dioxide 25 • Argon

Media should be defined when ordering.
Can be used for CA, N2, CO2, Ar. Not recommended for Helium (He). Do not use for Oxygen (O2)!

OUTLET

TECHNICAL GASES 2.0 & VACUUM/HIGH FLOW


Specifications
BROEN-LAB UniFlex™ outlets have an integrated media indication according to EN 13792.
Inlet/outlet: BROEN-LAB UniFlex™
Max. wall thickness: 20mm
Weight: 0.5 kg

Materials with media contact
Metals: Brass
Rubber: FKM
Plastics: PVDF
Lubricant: Perfluoropolyether based



Ordering information
For standard configuration please refer to the item numbers below.

- Other options:
- other media
 - flange
 - straight rear wall connection
 - wall thickness 20-40mm
 - double wall mounting
 - table mounting
 - additional support in case of a thin panel

For further information, please contact your BROEN-LAB representative.





With fixed nozzle



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
2.0 gases	2511 100 10XX-40	2511 110 10XX-40	2511 120 10XX-40	2511 130 10XX-40
Vacuum	2511 100 1028-40	2511 110 1028-40	2511 120 1028-40	2511 130 1028-40

With quick release coupling



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
2.0 gases	2511 102 10XX-40	2511 112 10XX-40	2511 122 10XX-40	2511 132 10XX-40
Vacuum	2511 102 1028-40	2511 112 1028-40	2511 122 1028-40	2511 132 1028-40

Replace **XX** with one of these media types:



21 • Compressed air



23 • Nitrogen



24 • Carbon dioxide



25 • Argon



26 • Helium



28 • Low vacuum (1000 to 1 mbar)

Delivered with 90° Rear Wall Connection (max wall thickness: 20mm)

FRONT CONTROL VALVE

TECHNICAL GASES 4.0 & OXYGEN

Specifications

Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792.
Inlet/outlet: BROEN-LAB UniFlex™
Weight: 0.7 kg
Test pressure: 1.5 x working pressure

Headwork: Needle
Opening/closing function: 3 x 360°
Max. working pressure: 16 bar / 232 psi

Headwork: Micro flow
Opening/closing function: 7.5 x 360°
Max. working pressure: 16 bar / 232 psi

Materials with media contact

Metals: Brass
Rubber: FKM
Plastics: PVDF
Lubricant: Perfluoropolyether based (Oxygen approved)


All components are ultrasonically cleaned for 99,99% clean gases

Ordering information

For standard configuration please refer to the item numbers below.


- Other options:
- a valve for front side mounting
 - a valve without mounting kit

For further information, please contact your BROEN-LAB representative.



Available with or without service valve according to product numbers listed below

Configuration options



Service valve

	Without service valve Needle headwork	Without service valve Micro flow headrwork	With service valve Needle headwork	With service valve Micro flow headrwork
	2511 527 10XX	2511 550 10XX	2511 528 10XX	2511 551 10XX

Replace **XX** with one of these media types:

- 21** • Compressed air
- 22** • Oxygen
- 23** • Nitrogen
- 24** • Carbon dioxide
- 25** • Argon
- 26** • Helium
- 27** • Dinitrogen monoxide

Delivered with mounting accessories

OUTLET

TECHNICAL GASES 4.0 & OXYGEN

Specifications
BROEN-LAB UniFlex™ outlets have an integrated media indication according to EN 13792.
Inlet/outlet: BROEN-LAB UniFlex™
Max. wall thickness: 20mm
Weight: 0.5 kg


Materials with media contact
Metals: Brass
Rubber: FKM
Plastics: PVDF
Lubricant: Perfluoropolyether based (Oxygen approved)

All components are ultrasonically cleaned for 99,99% clean gases

Ordering information
For standard configuration please refer to the item numbers below.

- Other options:
- other media
 - flange
 - removable nozzle
 - other standout
 - straight rear wall connection
 - wall thickness 20-40mm
 - table mounting
 - additional support in case of a thin panel

For further information, please contact your BROEN-LAB representative.



With fixed nozzle



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
	2511 108 10XX-40	2511 118 10XX-40	2511 128 10XX-40	2511 138 10XX-40

With quick release coupling



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
	2511 185 10XX-40	2511 186 10XX-40	2511 187 10XX-40	2511 190 10XX-40

Replace **XX** with one of these media types:

- 21** • Compressed air
- 22** • Oxygen
- 23** • Nitrogen
- 24** • Carbon dioxide
- 25** • Argon
- 26** • Helium
- 27** • Dinitrogen monoxide

Delivered with 90° Rear Wall Connection (max wall thickness: 20mm)

FRONT CONTROL VALVE

PURE 5.0 GASES

Specifications

Valve can be used for non-corrosive and inactive gases with purity coefficient up to 5.0 (99.999% pure gases).

Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792.
Inlet/outlet: BROEN-LAB UniFlex™
Leak rate: < 5x10⁻⁵ mbar x l/s [Helium]
Weight: 0.7 kg
Test pressure: 1.5 x working pressure

Headwork: Needle
Opening/closing function: 3 x 360°
Max. working pressure: 21 bar / 305 psi

Headwork: Micro flow
Opening/closing function: 7.5 x 360°
Max. working pressure: 16 bar / 232 psi

Shut-off Valve
Headwork: Stop
Opening/closing function: 1.5 x 360°
Max. working pressure: 21 bar / 305 psi

Materials with media contact

Metals: Brass
Rubber: FKM
Plastics: PVDF
Lubricant: Perfluoropolyether based (Oxygen approved)

All components are ultrasonically cleaned for 99,999% clean gases


Ordering information

For standard configuration please refer to the item numbers below.

Other options:

- a valve for front side mounting
- a valve without mounting kit
- inlet/outlet with female 1/4 NPT thread

For further information, please contact your BROEN-LAB representative.



Configuration options



	Shut-off valve	Needle valve	Micro flow valve
Fixed nozzle	2511 543 10XX	2511 542 10XX	2511 553 10XX

Replace **XX** with one of these media types:



21 • Compressed air



22 • Oxygen



23 • Nitrogen



24 • Carbon dioxide



25 • Argon



26 • Helium



27 • Dinitrogen monoxide

Delivered with mounting accessories

OUTLET
PURE 5.0 GASES

Specifications

Outlets can be used for non-corrosive and inactive gases with purity coefficient up to 5.0 (99,999% pure gases).

BROEN-LAB UniFlex™ outlets have an integrated media indication according to EN 13792.

Inlet: BROEN-LAB UniFlex™

Outlet: Female 1/4 NPT

Leak rate: < 5x10⁻⁵ mbar x l/s [Helium]

Max. wall thickness: 20mm

Weight: 0.5 kg


Materials with media contact

Metals: Brass

Rubber: FKM

Plastics: PVDF

Lubricant: Perfluoropolyether based (Oxygen approved)



Additional mounting information see page 30-33

All components are ultrasonically cleaned for 99,999% clean gases.

Ordering information



For standard configuration please refer to the item numbers below.

Other options:

- other media
- flange
- straight rear wall connection
- wall thickness 20-40mm
- table mounting
- additional support in case of a thin panel
- Swagelok® tube fittings
- inlet with female 1/4 NPT thread

For further information, please contact your BROEN-LAB representative.

Configuration options



	45°	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
	2511 109 10XX-40	2511 119 10XX-40	2511 129 10XX-40	2511 139 10XX-40

Replace **XX** with one of these media types:



21 • Compressed air



22 • Oxygen



23 • Nitrogen



24 • Carbon dioxide



25 • Argon



26 • Helium



27 • Dinitrogen monoxide

Delivered with 90° Rear Wall Connection (max wall thickness: 20mm)


FRONT CONTROL VALVE

PURE 6.0 GASES

Specifications

Valve can be used for non-corrosive and inactive gases with purity coefficient up to 6.0 (99.9999% pure gases).

Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792.
Headwork: Diaphragm type
Opening/closing function: 1.5 x 360°
Max. working pressure: 21 bar / 305 psi
Test pressure: 1.5 x working pressure
Inlet/outlet: Female 1/4" NTP thread
Leak rate: < 10-7 mbar x l/s [Helium]
Weight: 0.7 kg



Materials with media contact

Metals: Brass
Plastics: PVDF
Diaphragm: Stainless Steel Alloy

Ordering information

For standard configuration please refer to the item numbers below.

Media should be defined when ordering.

Other options:
- Swagelok® tube fittings

For further information, please contact your BROEN-LAB representative.

Configuration options



	Stop valve	Flow valve
All media	02 460.009-9023	02 462.009-9023

Media should be defined when ordering.
Delivered with mounting accessories

OUTLET

PURE 6.0 GASES

Specifications

BROEN-LAB CLEANLINE outlets for 6.0 gases can be used for mounting in the fume hoods in combination with BROEN-LAB CLEANLINE fittings with rear outlet connections.

Inlet connection: Female 1/4"NPT - male 1/2"NPT
(see ordering information for other options)

Outlet connection: Female 1/4"NPT
(see ordering information for other options)

Leak rate: < 10⁻⁷ mbar x l/s [He]

Weight: 0.5 kg

Materials with media contact

Metals: Brass or Stainless Steel



Swagelok® tube fittings
must be ordered separately

Configuration options



	45°	Straight	Bulk head union
Brass (powder coated, RAL 9023)	02869529-9023	02861529-9023	02566009-9023
Stainless Steel	02869533	02861533	

Delivered with mounting accessories.

FRONT CONTROL VALVE

CORROSIVE GAS

Specifications

Valve can be used for corrosive gases with purity coefficient up to 6.0 (99.9999% pure gases).

Handle: Metal, BROEN-LAB POLYCOAT powder coated with media indication according to EN13792.
Headwork: Diaphragm type
Opening/closing function: 1.5 x 360°
Max. working pressure: 16 bar / 232 psi
Test pressure: 1.5 x working pressure
Inlet/outlet: SS tubes for OD 8mm
Leak rate: < 10⁻⁷ mbar x l/s [Helium]
Weight: 0.7 kg



Materials with media contact

Metals: Stainless Steel (AISI 316L)
Plastics: PVDF
SS, AM350

Ordering information

For standard configuration please refer to the item numbers below.

Media should be defined when ordering.

Other options:
- a valve without mounting kit

For further information, please contact your BROEN-LAB representative.

To connect flexible pipes to this solutions, remember to offer three Swagelok® from page 29.
See how to install the solution on page 33.

Configuration options



	Without service valve
 Corrosive gases	2511 539 1099

Media should be defined when ordering.
Delivered with mounting accessories

OUTLET
CORROSIVE GAS

Specifications

Outlets can be used for corrosive gases with purity coefficient up to 6.0 (99.9999% pure gases).

BROEN-LAB UniFlex™ outlets have an integrated media indication according to EN 13792.

Inlet: OD 8mm / ID 6mm

Leak rate: < 10⁻⁷ mbar x l/s [Helium]

Max. wall thickness: 20mm

Weight: 0.6 kg

Materials with media contact

Metals: Stainless Steel (AISI 316L)

Plastics: PVDF

Ordering information


For standard configuration please refer to the item numbers below.

Media should be defined when ordering.

Other options:

- flange
- wall thickness 20-40mm
- additional support in case of a thin panel

For further information, please contact your BROEN-LAB representative.




Additional mounting information see page 30-33

To connect flexible pipes to this solutions, remember to offer three Swagelok® from page 29.

See how to install the solution on page 33.

With fixed nozzle









	90° S = 53 mm	90° S = 75 mm	90° S = 90 mm
 Corrosive gases	2511 150 10XX-40	2511 151 10XX-40	2511 152 10XX-40




Media should be defined when ordering.

Delivered with 90° Rear Wall Connection (max wall thickness: 20mm)

ACCESSORIES


Illustration	Description	Media	BROEN-LAB item no.
	Flange (wall thickness 0-20mm)	All media	1521 0579G
	Flange (wall thickness 20-40mm)	All media	1521 1789G
	Nut (wall thickness 0-20mm)	All media	1521 036
	Nut (wall thickness 20-40mm)	All media	1521 054
	Straight Rear Wall Connection	Water Burning gas Technical gases 2.0	1521 088
	90° Rear Wall Connection	Water Burning gas Technical gases 2.0	1521 052
	90° Rear Wall Connection	Special water	2511 492
	90° Rear Wall Connection	Technical gases 4.0 + 5.0	1521 0520
	Rear Wall support plate	All media	1521 092
	Test equipment for outlets	Technical gases 2.0	1521 177

ACCESSORIES


Illustration	Description	Media	BROEN-LAB item no.
	Rear wall connection for fume hoods with a double wall construction M = 49 to 61 mm M = 61 to 73 mm M = 73 to 85 mm M = 85 to 97 mm M = 97 to 109 mm	Water Burning gases Technical gases 2.0	2511 494 2511 497 2511 496 2511 495 2511 493
	Rear wall connection for fume hoods with a double wall construction M = 61 to 73 mm M = 73 to 85 mm M = 85 to 97 mm M = 97 to 109 mm	Special water	1530 479 1530 515 1530 534 1530 481
	Rear wall connection for fume hoods with a double wall construction M = 49 to 61 mm M = 61 to 73 mm M = 73 to 85 mm M = 85 to 97 mm	Technical gases 4.0 + 5.0	2511 4940 2511 4970 2511 4960 2511 4950
	Water jet filter pump for spouts (inlet: union nut G3/4, adapter to female G1/2 is included)	Water	1919 9009
	Water jet filter pump for straight outlets (inlet: union nut G3/4, adapter to female G3/8 is included)	Water	1530 4769
	Swagelok® Union OD10-OD8	Water Special water Burning gases Technical gases 2.0	96i0160
	Swagelok® Union OD10-OD8	Technical gases 4.0 + 5.0 + 6.0 Cleanline	96i01600

MOUNTING INSTRUCTIONS


FRONT CONTROL VALVES



1)




2)




MOUNTING INSTRUCTIONS

OUTLETS

Flange (optional)
Mandatory on soft wall materials



1)




2)




MOUNTING INSTRUCTIONS

MOUNTING OF OUTLETS INTO FUME HOODS WITH A DOUBLE WALL CONSTRUCTION


- 1) Insert the rear wall connection through the back wall



- 2) Fix the outlet to the back wall by tightening the nut



- 4) Position the outlet




- 5) Tighten Allen screws with Allen key, max 1Nm

PIPE INSTALLATION

WATER / SPECIAL WATER (SS MODEL) / BURNING GASES /
TECHNICAL GASES (2.0, 4.0 AND 5.0)

Valve

Outlet




BROEN-LAB UniFlex™ connection

Pipe installation can be done either with flexible hoses or with standard CU, PEX or SS pipes.

BROEN-LAB has a wide range of flexible pipes and hoses available with BROEN-LAB UniFlex™ connections as well as connection fittings to fit pipes of the following sizes:

- OD 8mm
- OD 10mm
- OD 12mm
- OD 15mm
- 3/8"

BROEN-LAB UniFlex™ flexible pipes and hoses



BROEN-LAB UniFlex™ connection fittings and adaptors
for CU, SS, PEX and PA pipes



PIPE INSTALLATION

SPECIAL WATER (PP MODELS ONLY)

Please note that PP valves and outlets for special water are delivered with push-in fittings.

BROEN-LAB recommends the following materials for pipe installation:

Valve

Outlet



Push-in fittings for OD10mm plastic tubing

BROEN-LAB recommends the following materials for pipe installation:

- PVDF (can be purchased at BROEN-LAB), LLDPE, PUR

PIPE INSTALLATION


CORROSIVE GASES

Please note that valves and outlets for corrosive gases are delivered in stainless steel and are prepared for orbital welding.

In case of outlet, welding process should take place before mounting into the fume hood.


Valve

Outlet



To connect flexible pipes to this solutions, remember to offer three Swagelok® from page 29.

See how to install the solution below.



NOTES

NOTES

HEADQUARTER

BROEN-LAB A/S

Drejervaenget 2
5610 Assens
Denmark

Tel. +45 6376 6376

lab@broen-lab.com
broen-lab.com

INTERNATIONAL OFFICES

BROEN-LAB Ltd.

Tel. +44 121 522 4515

lab@broen-lab.com
broen-lab.com

BROEN-LAB Singapore Pte. Ltd.

25 International Business Park
#04-60A German Centre
Singapore 609916

Tel. +65 6298 0662

lab@broen-lab.com
broen-lab.com

BROEN-LAB GmbH

Tel. +45 6376 6376

lab@broen-lab.com
broen-lab.com

BROEN-LAB Sverige AB

Tel. +45 6376 6376

lab@broen-lab.com
broen-lab.com

BROEN-LAB Inc.

15 Constitution Drive
Suite 122
NH 03110
Bedford
USA

Tel. +1 603 310 5089

lab@broen-lab.com
broen-lab.com

BROEN-LAB develops, manufactures and sells Laboratory Fittings, Emergency Shower Systems and Eye Wash Systems into a broad spectrum of laboratories and industrial locations; sectors include Pharmaceutical, Food&Beverage and Academia. Our expertise and product quality ensure optimal solutions compliant to all relevant international norms setting new standards in risk mitigation in modern research and development facilities.

BROEN-LAB offers solutions that ensure the functionality, compliance, hygiene, durability and safety are of the highest priority. For more than 50 years our products have been integrated into a wide variety of workplaces, laboratories, hospitals and industrial locations, with features that are assessed and designed to be flexible, durable and compliant offering a broad range of options tailored to each location. This inherent design integrity provides confidence in the solution adopted wherever in the world these are applied.

BROEN-LAB is a collaboration partner in the design and layout of your next laboratory and it's integrated safety.



BROEN-LAB A/S ISO 9001 certification

In September 1991 BROEN-LAB A/S was certified according to ISO 9001 as one of the first Danish companies. The certification was carried out by Bureau Veritas, Denmark, for our Danish site in Assens.

The quality management system of BROEN-LAB A/S now complies with detailed specifications laid down by the internationally acknowledged EN ISO 9001:2015. This certification will further contribute to reputation for quality and reliability of BROEN-LAB A/S.

BROEN-LAB A/S ISO 14001:2015 certification

In December 2022 we got certified according to ISO 14001 in Environmental certification. The new certification was carried out by Bureau Veritas, Denmark, for Danish site in Assens. The Environment ISO14001 management system of BROEN-LAB A/S is included in internationally acknowledged EN ISO 9001:2015. Both certifications furthermore contribute to reputation for quality and reliability of BROEN-LAB A/S.