


	Spare part number	Description
①	17481009	Outlet pipe with showerhead
②	17468009	Ballofix valve for ESS, short
③	17469009	Pull-rod with sign, Body and eye shower
④	17471009	3/4" pipe with fittings
⑤	17463009	Ballofix valve for eye shower
⑥	17466009	Screws and O-rings
⑦	17464009	Push arm
⑧	17462009	Bowl for eyeshower
⑨	17461009	Eye shower unit for ESS
⑩	17478009	Floor bracket for body & eye shower

EN-Label: 96M0375



Description: Body & eye shower with bowl; Freestanding	
Media: Water	
Material:	
Specification: Premium	
Comments:	
Drawing type: Assembly drawing	Weight: 11.844 Kg
All dimensions made with  are inspection dimensions	
Unspecified linear tolerances:	Unspecified angular tolerances:
Designed by: HEM	Drawn by: TLA
Released on TM No. 6979	Date: 09-08-17
Format: A2	Drawing No.: 17556009SP

Rev.	Pos.	Date	TM No.	Design by	Drawn by	Change:

**BROEN-LAB** This drawing is our property and cannot be given to a third party, copied or otherwise used without our permission. BROEN-LAB A/S

