

Area of Application:  
Eye and body shower  
Water  
Min. working pressure 2 bar  
Adjustable body shower head



Min. working pressure 1,5 bar

Materials in media contact:  
Brass  
Rubber, EPDM  
PP  
Silicone based lubricant



Inlet / outlet connection:  
Inlet G1" female  
Water outlet DIA 1 1/2" (Ø40mm) drain pipe

Cabinet materials:  
AISI 316L Stainless steel  
Brushed.  
Drain pan: 18 gauge Stainless steel  
Cabinet: 16 gauge Stainless steel


Signs & instructions:  
(All signs to be mounted on site)



Eye wash compartment



Heavy duty stainless steel handle with black vinyl grip to activate body shower



Description: Eye and body shower cabinet

						Media: Water	
						Comments:	
						Drawing type: Customer drawing	Weight: 19.995 Kg
Rev.	Pos.	Date	TM No.	Design by	Drawn by:	Designed by: NSC	Released on TM No.
						Drawn by: CGL	Date (DD-MM-YY): 1-02-24
						Format: A3	Drawing No.: 17157009C

**BROEN-LAB**

This drawing is our property and cannot be given to a third party, copied or otherwise used without our permission. BROEN-LAB A/S.