

4


3

2

1


	Spare part number	Description
①	17479009	Body shower outlet, ceiling mounted
②	17468009	Ballofix valve for ESS, short ⑦
③	17488009	Pull rod incl. body shower sign
④	17480009	Ceiling flange connection
⑤	19138000	Mounting kit
⑥	17466009	Screws and O-rings
⑦	17001910	Body shower sign

EN-Label: 96M0373



Description:	Body shower, ceiling mounted
Media:	Water
Material:	
Specification:	Premium

Comments:	
Drawing type:	Assembly drawing
Weight:	1.540 Kg

All dimensions made with  are inspection dimensions  
 Unspecified linear tolerances:      Unspecified angular tolerances:

Rev.	Pos.	Date	TM No.	Design by	Drawn by	Change:

Designed by:	Drawn by:	Released on TM No.	Date:
HEM	SJE	6979	09-08-17

**BROEN-LAB** This drawing is our property and cannot be given to a third party, copied or otherwise used without our permission. BROEN-LAB A/S

Format: Drawing No.: **17580009SP**

4

3

2

1