

Area of Application:
Turret
Vacuum & Highflow

Max. working pressure 16 bar

Details:
High flow/Vacuum headwork open/close 1,5x360°
3-Zone colour, acc. to EN13792:2002





Fixed nozzle

Materials in media contact:
Brass
PP
Rubber, FKM
Perfluoropolyether based lubricant
SPX (Soft PEX)
PVDF

Mounting:
M30x1,5 nut with locking pins, O-ring and washer (See 94G0232)

Inlet / outlet connection:
DN8 hose with G1/2" female connection
DIN 12 898 hose nozzle

Polypropylene handle



MEDIA CODE XX FOR VACUUM *Model shown for V		
Media	Media description	Media code XX
V	Vacuum Low	28
VF	Vacuum Fine	29
VH	Vacuum High	30

Description: 1-Valve, Turret, 1-Hole, Fix. Nozzle, M30x1,5mm, Hose

								Media:	Vacuum & Highflow		
								Comments:			
								Drawing type:	Customer drawing	Weight:	1.014 Kg
Rev.	Pos.	Date	TM No.	Design by	Drawn by:	Change:		Designed by:	Drawn by:	Released on TM No.	Date (DD-MM-YY):
								VJE	MWA		10-08-20
BROEN-LAB This drawing is our property and cannot be given to a third party, copied or otherwise used without our permission. BROEN-LAB A/S.							Format:	Drawing No.: 253200210XX-52E			
							A3				