

Area of Application:
Table/Wall mounted
Water

Max. working pressure 10 bar
Min. temperature 5°C
Max. temperature 90°C



Details:
Ceramic headwork open/close 180°
3-Zone colour, acc. to EN13792:2002
Swivel spout
1/2" removable nozzle

Materials in media contact:
Brass
Rubber, EPDM
Silicone based lubricant
SPX (Soft PEX)
POM
Ceramic


Mounting:
1/2" nut with locking pins o-ring and washer (See 94G0232)

Inlet / outlet connection:
DN8 hose with G3/8 female inlet
DIN 12 898 hose nozzle

Metal handle



WALL OR BENCH MOUNTED



Description: 2-Hand.Mix, Table/Wall mounted, 2-Hole, Rem. Nozzle, M30x1,5mm, Hose

							Media:	Water		
							Comments:	150mm center distance		
							Drawing type:	Customer drawing	Weight:	3.098 Kg
Rev.	Pos.	Date	TM No.	Design by	Drawn by	Change:	Designed by:	Drawn by:	Released on TM No.	Date:
							BJP	MWA		24-05-2019
BROEN-LAB							Format:	Drawing No.:	253606110XX-52D	
							A3			

This drawing is our property and cannot be given to a third party, copied or otherwise used without our permission. BROEN-LAB A/S.

MEDIA CODE XX FOR WATER
TWO HANDLE MIXER

Media	Media description	Media code XX
WPC WPH	Water potable cold / hot	39
WNC WCH	Water non-potable cold / hot	40