



Area of Application:
Body Shower & Eye Wash, Wall mounted
Min. operation pressure: 1,5 Bar
Max. operation pressure: 10 Bar
Min. flow 60 l/min
Product temperature range 1°C-90°C
Operating water temperature range 15°C-37°C (ideal 20°C-25°C)

Details:
BROEN-LAB Stainless steel Glass-blasted surface
Signs: Safety sign package Combi Shower (17001920)
Drain: Male ISO 228/1 - G5/4"
Design acc. EN15154-1 & EN15154-5 (Class 2 - 60 to 100 l/min)

Materials in media contact:
Main: SS316L, PTFE, Rubber EPDM, Silicone based lubricant, Foodlube MultiPasta
Eye Shower: PP, PA6, PPH, TPE, Rubber EPDM, Silicone based lubricant
Body Shower head: PPO, PP, Rubber EPDM, Silicone based lubricant
Flow reduction unit: POM, EPDM

Mounting:
Wall bracket Ø5,5 x 6 (60°) (See 94G0298)

Inlet / outlet connection:
Female ISO 228/1 - G1"
Body & Eye shower Spray Heads



*Optional inlet - ** Blingplug for unused Inlet

Description: Combination Shower, Wall, 316L, Glass Blasted, Glass Blasted									
							Media: Water		
							Comments:		
							Drawing type: Customer drawing		Weight: 14.606 Kg
Rev.	Pos.	Date	TM No.	Design by	Drawn by	Change:	Designed by: BDA	Drawn by: ACH	Released on TM No. 10602
BROEN-LAB							Date (DD-MM-YY): 23-06-26		
							Format: A3	Drawing No.: 28433031021000	

This drawing is our property and cannot be given to a third party, copied or otherwise used without our permission. BROEN-LAB A/S.