



Area of Application:
 Body Shower & Eye Wash, Wall mounted
 Min. operation pressure: 1,5 Bar
 Max. operation pressure: 10 Bar
 Min. flow 60 l/min
 Product temperature range 1°C-90°C
 Operating water temperature range 15°C-37°C (ideal 20°C-25°C)

Details:
 BROEN-LAB Stainless steel Glass-blasted surface
 Signs: Safety sign package Combi Shower (17001920)
 Drain: Male ISO 228/1 - G5/4"
 Design acc. EN15154-1 & EN15154-5 (Class 2 - 60 to 100 l/min)

Materials in media contact:
 Main: SS316L, PTFE, Rubber EPDM, Silicone based lubricant, Foodlube MultiPasta
 Eye Shower: PP, PA6, PPH, TPE, Rubber EPDM, Silicone based lubricant
 Body Shower head: PPO, PP, Rubber EPDM, Silicone based lubricant
 Flow reduction unit: POM, EPDM

Mounting:
 Wall bracket Ø5,5 x 6 (60°) (See 94G0298)

Inlet / outlet connection:
 Female ISO 228/1 - G1"
 Body & Eye shower Spray Heads



*Optional inlet - ** Blingplug for unused Inlet

Description: Combination Shower, Wall, 316L							Media: Water				
							Comments:				
							Drawing type: Customer drawing				
							Weight: 14.650 Kg				
Rev.	Pos.	Date	TM No.	Design by	Drawn by	Change:	Designed by: MJO	Drawn by: MJO	Released on TM No.	Date (DD-MM-YY): 27-08-25	
BROEN-LAB							This drawing is our property and cannot be given to a third party, copied or otherwise used without our permission. BROEN-LAB A/S.		Format: A3		Drawing No.: 28433031021006